

Table 1

# Cells

Cell	Receptor
5.25	CXCR4, CD4, CCR5 (not expressed well) BONZO
5.25.Luc4.M7	CD4, CCR5, BONZO
HOS.CD4.CCR5	CD4, CCR5
HOS.CD4.CXCR4	CD4, CXCR4
HOS.CD4	CD4, low level expression of CCR5 and CXCR4
HOS HT4 R5 GFP wt	CD4, CXCR4, CCR5
HOS.CD4.CCR5.GFP.M7#6*	CD4, CXCR4, CCR5
P4.CCR5	CD4, CXCR4, CCR5
U87.CD4	CD4
U87.CD4 R5	CD4, CCR5
U87.CD4 X4	CD4, CXCR4
MT2	CD4, CXCR4
MT4	CD4, CXCR4
PM1	CD4, CXCR4, CCR5
CEM NKr CCR5	CD4, CXCR4, CCR5

**Table 2** Representative viruses and reagents

Viruses	Envelope <sup>a</sup>	Source
89.6, SF2	R5-X4/SI/B	ARRRP <sup>b</sup>
92BR014, 92US076	R5-X4/SI/B	ARRRP
JR-CSF, 91US005	R5/NSI/B	ARRRP
91US054	SI/B	ARRRP
NL43, MN, ELI	X4/B	ARRRP
92HT599	X4	ARRRP
92UG031	R5/NSI/A	ARRRP (IN-HOUSE)
92TH014, 92TH026	R5/NSI/B	ARRRP (IN-HOUSE)
92BR025, 93MW959	R5/SI/C	ARRRP (IN-HOUSE)
92UG035	R5/NSI/D	ARRRP (IN-HOUSE)
92TH022, 92TH023	R5/NSI/E	ARRRP (IN-HOUSE)
93BR020	R5-X4/SI/F	ARRRP (IN-HOUSE)
Antibodies	Epitope	SOURCE
Mabs 2F5, 1577	gp41 TM	ARRRP
Mabs IG1b12, 2G12, 17b, 48D	gp120 SU	ARRRP
Neutralization sera #2, HIV-IG	Polyclonal	ARRRP
Entry inhibitors	Target	Source
CD4-IG	gp120 SU	Genentech
CD4-IGG2	gp120 SU	Adarc
SCD4	Sigma	Progenics
T20 (DP178)	gp41 TM	Trimeris
Rantes, MIPl a/b	CCR5	SIGMA/ARRRP
SDF1a/b	CXCR4	SIGMA/ARRRP
AMD 3100	CXCR4	AnorMed
Dextran sulfate, Heparin	Non-specific	Sigma

<sup>a</sup>R5 (CCR5 co-receptor), X4 (CXCR4 co-receptor)

SI (syncytium inducing), NSI (non-syncytium inducing), A,B,C,D,E,F (envelope clade designation)

<sup>b</sup>AIDS Research and Reference Reagent Program

# Table 3

## Primers Tested for the Amplification of HIV Envelope

### RT PRIMERS

RT env\_N3 5'-GGA GCA TTT ACA AGC AGC AAC ACA GC-3'  
 RT env 9720 5'-TTC CAG TCA VAC CTC AGG TAC-3'  
 RT env 9740 5'-AGA CCA ATG ACT TAY AAG G-3'

### 5' PCR PRIMERS

5'env 5'-GGG CTC GAG ACC GGT CAG TGG CAA TGA GAG TGA AG-3'  
 5'env\_Xho/Pin 5'-GGG CTC GAG ACC GGT GAG CAG AA-3 ACA GTG GCA ATG A-3'  
 5'env\_START 5'-GGG CTC GAG ACC GGT GAG CAG AA-3 ACA GTG GCA ATG -3'

### 3' PCR PRIMERS

3' env 5'-GGG TCT AGA ACG CGT TGC CAC CCA TCT TAT AGC AA-3'  
 3'env\_Xba/Mlu 5'-GGG TCT AGA ACG CGT CCA CTT GCC ACC CAT BTT ATA GC-3'  
 3'env\_STOP 5'-GGG TCT AGA ACG CGT CCA CTT GCC ACC CAT BTT A-3'  
 3' delta CT 5'-GAT GGT CTA AGA CGC TGT TCA ATA TCC CTG CCT AAC TC-3'

All Experiments are located in Virologic Book number 0188

# Table 4 (Panel 1)

Anti-HIV Drugs			
Drug/Compound	Generic Name	Trademark	Manufacturer
<b>RT Inhibitors (NRTI, nucleotide analogs)</b>			
AZT, ZDV	Zidovudine	Retrovir	Glaxo/Wellcome
3TC	Lamivudine	Epivir	Glaxo/Wellcome
AZT + 3TC		Combivir	Glaxo/Wellcome
d4T	Stavudine	Zerit	Bristol-Myers/Squibb
ddl	Didanosine	Videx	Bristol-Myers/Squibb
ddC	Zalcitabine	Hivid	Hoffman La Roche
1592U89	Abacavir	Ziagen	Glaxo/Wellcome
AZT + 3TC + 1592U89		Trizivir	Glaxo/Wellcome
(-)FTC (5-fluoro-3TC; Corviracil)	Emtricitabine		Triangle Pharmaceuticals
(-)FTC + (+)FTC (50:50)	Racimir		QuadPharma
DAPD (DXG active)	Amdoxovir		Triangle Pharmaceuticals
F-ddA (2-fluoro-ddA)	Lodenosine		MedImmune Oncology (US Bioscience)
BCH-10652, dOTC (2-deoxy-3-oxa-4-thiocytidine)			BioChem Pharma, Inc.
D-d4FC			Triangle Pharmaceuticals (Schinazi)
<b>RT Inhibitors (NNRTI, nucleotide analogs)</b>			
bis-POC PMPA (GS-4331)	Tenofovir		Gilead Sciences
bis-POM PMEa (GS-840)	Adefovir dipivoxil		Gilead Sciences
<b>RT Inhibitors (NNRTI, non-nucleosides)</b>			
BI-RG-587	Nevirapine	Viramune	Boehringer/Ingelheim (Roxanne)
BHAP PNU-90152T	Delavirdine	Rescriptor	Pharmacia & Upjohn
DMP 266 (L-743,726)	Efavirenz	Sustiva	Dupont Pharmaceuticals (Avid)
MKC442 (Coactinon)	Emivirine		Triangle/Mitsubishi Kasei
AG-1549 (S1153) (on hold)	Capravirine		Agouron Pharmaceuticals
PNU-142721			Pharmacia & Upjohn
DPC-961, -963, -083, -08?			DuPont Pharmaceuticals
SJ-3366	Also entry inhibitor?		Samjin Pharmaceuticals
BHAP PNU-87201	Ateviridine		Upjohn
GW420867X (quinoxaline)	(2 <sup>nd</sup> gen. HBY 097)		Glaxo/Wellcome (Hoechst Bayer)
TMC 120 (R147681)			Tibotec
TMC 125 (R165335)			Tibotec
R86183	tivirapine		Janssen Pharmaceuticals
Calanolide A			Sarawak Medichem Pharmaceuticals
<b>Protease Inhibitors (PRI)</b>			
Ro 31-8959	Saquinavir-(hgc) Saquinavir-(sgc)	Invirase Fortivase	Hoffman-La Roche
MK-639 (L-735,524)	Indinavir	Crixivar	Merck Research Laboratories
ABT-538 (A-84538)	Ritonavir	Norvir	Abbott Laboratories
AG1343	Nelfinavir	Viracept	Agouron Pharmaceuticals
141W94 (VX-478)	Ampranavir	Agenerase	Glaxo-Wellcome/Vertex
ABT-378/r	Lopinavir/ritonavir	Kaletra	Abbott Laboratories
BMS 232,632 (aza-peptide)			Bristol-Myers-Squibb
PNU-140690	Tipranavir		Pharmacia & Upjohn
DMP 450 (cyclic urea)	Mozenavir		Triangle/Avid (ph I/II)
TMC 126 (Erickson's compound)			Tibotec
G/W433908 (VX-175)	amprenavir pro-drug		Glaxo/Wellcome/Vertex
L756,423 (on hold)			Merck
PD-178390 (dihydropyrone)			Parke Davis (Boehringer-Ingelheim)
? new candidate			Roche
DPC 681 and 684			DuPont Pharmaceuticals
AG-1776 (JE-2147 = KNI-764)			Agouron Pharmaceuticals
<b>Envelope/Receptor Inhibitors</b>			
T-20 (gp41)	Pentafuside		Trimeris Pharmaceuticals
T-1249 (gp41)			Trimeris Pharmaceuticals
D-peptide inhibitor (gp41) small mol.	SCH-C		Schering-Plough
AMD-3100 (CXCR4)	(bicyclam)		AnorMED
AMD-8664 (CXCR4)	(macrocyclam)		AnorMED
ALX40-4C (CXCR4)			U. PA
FP21399			Fuji Pharmaceuticals
PRO 542 (gp120)	CD4lgG2		Progenics Pharmaceuticals
PRO-140 (CCR5)	MAb CCR5		Progenics Pharmaceuticals
T-22 (CXCR4)	(peptide, 18-mer)		
Met-SDF-1 (CXCR4)			
TAK 779 (CCR5 antagonist)			Takeda
AOP-Rantes (CCR5)			Gryphon Sciences

# Table 4 (Panel 2)

Rantes 9-68 (CCR5)			
CCR5 antagonists	4-(piperidin-1-yl) butane class		Merck
$\alpha$ -Immunokine-NNS03 (CCR5, CXCR4)	$\alpha$ -cobratoxin		PhyloMed Corp.
Integrase Inhibitors			
AR-177	Zintevir		Aronex Pharmaceuticals
Diketo acids			Merck Research Laboratories
Nucleocapsid Inhibitors			
RB 2121	cyclic peptide p7 mimic		(see PNAS 96:4886-4891 (1999))
CI-1012			Achelion Pharmaceuticals
RNase H Inhibitor			
SP1093V (BBNH Fe+3 derivative)			(Parniak)

FDA approved drugs are shown in boldface, red = discontinued development, blue = not sure about development status

# Table 4 (Panel 3)

Generic Name (abbreviation)	Brand Name	Firm	FDA Approval Date
zidovudine, AZT	Retrovir	Glaxo Wellcome	March 87
didanosine, ddl	Videx	Bristol Myers-Squibb	October 91
zalcitabine, ddC	Hivid	Hoffman-La Roche	June 92
stavudine, d4T	Zerit	Bristol Myers-Squibb	June 94
lamivudine, 3TC	Epivir	Glaxo Wellcome	November 95
saquinavir, SQV, hgc	Invirase	Hoffman-La Roche	December 95
saquinavir, SQV, sgc	Fortovase	Hoffman-La Roche	November 97
ritonavir, RTV	Norvir	Abbott Laboratories	March 96
indinavir, IDV	Crixivan	Merck & Co., Inc.	March 96
nevirapine, NVP	Viramune	Boehringer Ingelheim	June 96
nelfinavir, NFV	Viracept	Agouron Pharmaceuticals	March 97
delavirdine, DLV	Rescriptor	Pharmacia & Upjohn	April 97
ZDV+3TC	Combivir	Glaxo Wellcome	September 97
efavirenz, EFV	Sustiva	DuPont Pharmaceuticals	September 98
abacavir, ABC	Ziagen	Glaxo Wellcome	February 99
amprenavir	Agenerase	Glaxo Wellcome	April 99
lopinavir/ritonavir	Kaletra	Abbott	September 2000
ZDV+3TC+ABC	Trizivir	GlaxoSmithKline	November 2000

5.

drug 1: L83  
drug 2: AMD

Year	1970	1971	1972	1973	1974	1975	1976	1977	1978	1979	1980	1981	1982	1983	1984	1985	1986	1987	1988	1989	1990	1991	1992	1993	1994	1995	1996	1997	1998	1999	2000	2001	2002	2003	2004	2005	2006	2007	2008	2009	2010	2011	2012	2013	2014	2015	2016	2017	2018	2019	2020	2021	2022	2023	2024	2025	2026	2027	2028	2029	2030	2031	2032	2033	2034	2035	2036	2037	2038	2039	2040	2041	2042	2043	2044	2045	2046	2047	2048	2049	2050	2051	2052	2053	2054	2055	2056	2057	2058	2059	2060	2061	2062	2063	2064	2065	2066	2067	2068	2069	2070	2071	2072	2073	2074	2075	2076	2077	2078	2079	2080	2081	2082	2083	2084	2085	2086	2087	2088	2089	2090	2091	2092	2093	2094	2095	2096	2097	2098	2099	2100
1970	1971	1972	1973	1974	1975	1976	1977	1978	1979	1980	1981	1982	1983	1984	1985	1986	1987	1988	1989	1990	1991	1992	1993	1994	1995	1996	1997	1998	1999	2000	2001	2002	2003	2004	2005	2006	2007	2008	2009	2010	2011	2012	2013	2014	2015	2016	2017	2018	2019	2020	2021	2022	2023	2024	2025	2026	2027	2028	2029	2030	2031	2032	2033	2034	2035	2036	2037	2038	2039	2040	2041	2042	2043	2044	2045	2046	2047	2048	2049	2050	2051	2052	2053	2054	2055	2056	2057	2058	2059	2060	2061	2062	2063	2064	2065	2066	2067	2068	2069	2070	2071	2072	2073	2074	2075	2076	2077	2078	2079	2080	2081	2082	2083	2084	2085	2086	2087	2088	2089	2090	2091	2092	2093	2094	2095	2096	2097	2098	2099	2100	

	7			8			9			10			11			12		
No Drug RLU L83 RLU AMD RLU	R5 30	X4 1.217 814 42	R5:X4 0.02	R5 49	X4 1,128 934 37	R5:X4 0.04	R5 42,206 90 33,192	X4 70	R5:X4 602 94	R5 45	X4 1,159 761 33	R5:X4 0 04	R5 44,589 60 30,510	X4 97	R5:X4 459 68	R5 32	X4 29,118 22,279 31	R5:X4 0.00
No Drug RLU L83 RLU AMD RLU	R5 100	X4 5,204 3,058 34	R5:X4 0 02	R5 11,299 22 8,232	X4 52	R5:X4 217.29	R5 1,273 26 1,129	X4 20	R5:X4 63 65	R5 7,375 10,856 3,961	X4 4	R5:X4 0 68	R5 4,397 7,130 42 1,656	X4 7,130 3,874 1,285	R5:X4 0 62	R5 16,115 68 11,980	X4 50	R5:X4 322.30
No Drug RLU L83 RLU AMD RLU	R5 53	X4 974 828 95	R5:X4 0.05	R5 177 32 361	X4 32	R5:X4 5.53	R5 24,739 31 21,767	X4 33	R5:X4 749.67	R5 612 48 479	X4 45	R5:X4 13 60	R5 2,997 9,695 22 925	X4 9,695 7,261 35	R5:X4 0 31	R5 416 28 597	X4 31	R5:X4 13.42
No Drug RLU L83 RLU AMD RLU	R5 45	X4 338 207 29	R5:X4 0.13	R5 714 31 366	X4 45	R5:X4 15 87	R5 11,029 18 8,390	X4 18	R5:X4 612 72	R5 2,908 31 2,306	X4 43	R5:X4 67 63	R5 13,997 73 12,139	X4 52	R5:X4 269.17	R5 24,377 59 16,045	X4 34	R5:X4 716.97
No Drug RLU L83 RLU AMD RLU	R5 9,254 43 2,219	X4 23,846 9,150 40	R5:X4 0 39	R5 5,249 38 764	X4 20,393 8,290 38	R5:X4 0 26	R5 470 37 310	X4 39	R5:X4 12 05	R5 19,175 96 19,638	X4 126 183 35	R5:X4 152 18	R5 513 35 500	X4 31	R5:X4 16.55	R5 3,264 30 3,038	X4 28	R5:X4 116 57
No Drug RLU L83 RLU AMD RLU	R5 3,853 39 2,108	X4 77	R5:X4 50 04	R5 411 22 253	X4 34	R5:X4 12 09	R5 7,857 42 4,595	X4 41	R5:X4 191 63	R5 25,437 35 14,749	X4 38	R5:X4 669.39	R5 17,443 40 13,708	X4 45	R5:X4 387 62	R5 16,707 104 12,486	X4 38	R5:X4 439 66
No Drug RLU L83 RLU AMD RLU	R5 263 33 143	X4 38	R5:X4 6 92	R5 3,890 44 986	X4 9,577 6,140 45	R5:X4 0 41	R5 2,089 33 1,105	X4 1,698 1,224 160	R5:X4 1 23	R5 475 30 159	X4 43	R5:X4 11.05	R5 8,475 30 6,862	X4 38	R5:X4 223 03	R5 4,107 28 2,571	X4 33	R5:X4 124.45
No Drug RLU L83 RLU AMD RLU	R5 443 33 451	X4 13,452 12,051 51	R5:X4 0 03	R5 34 128,238 85,018 45	X4 R5:X4 0.00	R5 22	X4 62	R5:X4 R5:X4	R5 45	X4 33	R5:X4 R5:X4	R5 43	X4 42	R5:X4 R5:X4	R5 46	X4 35	R5:X4 R5:X4	

Table 4B

	3		
	R5 cells	X4 cells	R5:X4
no drug			
R5 inhibitor			
X4 inhibitor			
%inhib by R5 inhibitor			
%inhib by X4 inhibitor			

	21		
	R5	X4	R5:X4
no drug			
R5 inhibitor			
X4 inhibitor			
%inhib by R5 inhibitor			
%inhib by X4 inhibitor			

[illegible]

**DUAL or MIXED  
DEAD**

	76		
	R5	X4	R5:X4
no drug	14,982	12,020	1
R5 inhibitor	111	10,839	
X4 inhibitor	8,580	3,384	
%inhb by R5 inhibitor	99	10	
%inhb by X4 inhibitor	43	72	

		95	
	R5	X4	R5:X4
no drug	43	42	
R5 inhibitor			
X4 inhibitor			
%inhib by R5 inhibitor			
%inhib by X4 inhibitor			



Table 4C

Plate Repeat End time Start temp End temp. BarCode  
1 1 6:39:36 PM 21.6 21.7 N/A

0.5 CPS (CPS)  
18

26648	54	6970	37406	7158	38	1386	930	100	1184	112	30488
19248	262	72	18972	3478	3946	5010	46	14	11004	7474	60
106	9038	4002	32	238	26	976	34	38	48	9736	32
38	2000	528	16	16	36	298	52	14	54	52	38
122	2984	7264	40	38	20	23344	15340	44	156	40	32
36	52	28	40	28	40	62	32	38	22	44	28
22	32	28	9858	68	3802	46	11470	1958	42	42	24
30	64	36	3846	3390	57858	12620	126186	68	34	38	36

Plate Repeat End time Start temp End temp. BarCode  
2 1 6:41:51 PM 21.5 21.6 N/A

0.5 CPS (CPS)  
34

26294	56	6728	38882	7514	24	1048	1326	40	1134	82	27748
18956	190	72	20458	4418	4714	5398	58	26	10708	6786	40
30	15002	4766	42	246	44	972	30	28	42	9654	30
40	2620	808	32	28	86	378	38	22	32	52	30
92	5434	9638	36	42	38	24348	25446	34	96	22	24
44	44	14	46	46	62	92	36	44	54	46	48
30	44	28	8346	142	3468	30	7684	1438	44	34	42
36	54	30	3784	3840	65330	14284	130290	56	32	46	34

Plate Repeat End time Start temp End temp. BarCode  
3 1 6:44:06 PM 21.6 21.6 N/A

0.5 CPS (CPS)  
16

18590	18	4306	23902	5386	30	924	894	32	660	48	23382
8698	148	34	14088	2880	3142	3160	40	46	7616	2842	32
30	10252	3542	38	172	34	842	30	32	44	7616	34
28	2396	370	28	26	66	172	32	22	28	36	26
64	3784	4822	32	28	24	6188	5702	34	120	34	24
16	28	62	26	28	30	38	38	38	44	62	34
38	52	40	7020	68	1798	62	7324	1076	50	42	32
32	50	70	3824	3138	53670	13088	104608	32	38	50	32

Plate Repeat End time Start temp End temp. BarCode  
4 1 6:46:20 PM 21.7 21.8 N/A

0.5 CPS (CPS)  
24

19926	46	5086	27182	5550	26	704	974	32	862	44	21176
11214	158	42	15882	3026	2694	2956	26	20	5294	4906	44
56	11426	3252	82	270	52	814	36	40	64	6906	32
32	2916	556	42	30	64	242	38	24	36	42	40
72	3220	7618	36	38	44	12112	10878	32	246	32	42
32	48	22	30	28	24	136	30	56	34	24	24
32	38	40	8844	86	2386	54	4956	1372	50	46	28
22	34	40	3498	2754	59808	11014	65428	34	40	24	18

Plate Repeat End time Start temp End temp. BarCode  
5 1 6:48:35 PM 21.9 21.9 N/A

0.5 CPS (CPS)  
38

38	28	40	30	38	234	46	30	46	30	34	22
26	20	32	976	190	58	46	24	28	2298	1290	26
42	3260	420	50	36	34	102	26	34	24	40	24
28	378	88	32	34	32	40	42	54	26	30	48
36	62	11690	38	44	32	42	42	26	36	46	42
46	22	44	38	58	32	38	30	32	36	38	12
22	38	52	40	42	26	34	50	200	36	34	36
36	80	52	40	36	40	52	54	22	40	42	36

Plate Repeat End time Start temp End temp. BarCode  
6 1 6:50:50 PM 21.8 21.9 N/A

0.5 CPS (CPS)  
30

26	28	40	34	40	50	38	44	34	36	40	40
42	28	32	1068	114	38	22	46	30	2268	1280	40
36	3508	656	48	40	34	88	30	34	56	30	36
28	166	56	28	18	40	18	38	32	34	26	32
32	114	5732	32	42	34	38	34	50	34	38	52
36	22	36	40	30	20	30	26	42	40	26	36
46	24	26	32	60	28	32	40	120	36	20	36
32	58	58	32	24	40	50	36	12	38	40	18

Talk to D

Plate Repeat End time Start temp End temp BarCode  
1 1 ##### 21 7 21 8 N/A

0.5 CPS (CPS)  
32

24	14218	56	21586	28	4034	36	52	42088	40	43332	38
4780	174	71276	730	48	326	28	12022	1288	7198	4314	17856
69452	15306	4008	44	792	180	50	182	25292	718	3572	424
72	3922	546	128	44	66	30	900	11984	3194	16106	23794
30	88	18820	40438	4882	105946	8466	4934	470	20386	548	3376
3420	934	344	6268	5012	170	3546	504	8164	22214	15146	19592
344	82	5344	5330	6710	1880	338	4286	2112	466	7752	4078
580	1018	516	318	302	40	414	32	20	52	48	42

Plate Repeat End time Start temp End temp BarCode  
2 1 ##### 21 9 21 9 N/A

0.5 CPS (CPS)  
22

44	14062	36	21726	32	3942	24	46	42324	50	45846	26
4324	124	63502	492	44	242	172	10576	1258	7552	4480	14374
66204	14658	3568	56	524	282	56	172	24186	506	2422	408
128	3526	354	168	20	54	60	528	10074	2622	11888	24960
40	78	17378	38076	5944	89520	10042	5564	470	17964	478	3152
4118	638	336	4110	3624	222	4160	318	7550	28660	19740	13822
438	114	3554	3384	5470	1892	188	3494	2066	484	9198	4136
440	940	466	282	292	38	472	36	24	38	38	50

Plate Repeat End time Start temp End temp BarCode  
3 1 ##### 21.9 22 N/A

0.5 CPS (CPS)  
36

36	38	28	28	38	28	36	42	114	40	60	46
32	34	80	26	16	36	30	18	26	30	48	80
126	106	34	42	30	46	22	26	32	58	8	28
26	46	42	38	34	20	18	34	24	36	100	60
28	146	15756	76	32	26	42	52	36	106	34	24
34	50	20	24	32	18	38	18	36	38	50	58
38	38	28	36	80	50	38	44	34	38	30	28
40	30	38	30	24	22	38	36	34	36	22	26

Plate Repeat End time Start temp End temp BarCode  
4 1 ##### 21.9 22 N/A

0.5 CPS (CPS)  
16

24	70	30	26	26	40	18	32	66	30	60	34
24	46	118	34	26	40	34	26	26	50	36	56
132	116	84	34	34	48	34	38	30	38	36	28
32	30	32	28	26	40	46	28	12	26	46	58
34	84	12036	94	38	34	44	24	38	86	36	36
24	28	42	44	24	28	40	26	48	32	30	150
62	24	58	32	66	26	28	44	32	22	30	28
44	80	60	28	40	28	28	34	32	36	18	28

Plate Repeat End time Start temp End temp BarCode  
5 1 ##### 22 22 N/A

0.5 CPS (CPS)  
34

38	11132	26	6696	34	3960	44	40	29548	40	30628	34
762	48	39888	88	40	30	36	8416	1262	4096	1598	12372
52058	8050	2524	60	450	168	32	340	22410	534	1046	512
106	2080	172	170	38	38	30	444	7478	2458	10024	15470
30	46	16428	25792	4240	45094	2092	630	334	17130	388	2730
2498	732	358	5044	3236	202	2124	292	4806	13736	13012	12386
290	30	3986	886	4584	240	110	940	1286	162	6768	2204
456	816	438	294	260	44	524	40	34	22	38	20

Plate Repeat End time Start temp End temp BarCode  
6 1 ##### 22 22.1 N/A

0.5 CPS (CPS)  
32

40	13240	38	7756	38	3368	24	40	36836	48	30392	36
914	34	44702	64	50	20	24	8048	996	3826	1714	11588
50234	9110	1574	46	536	270	20	382	21124	424	804	682
200	1888	240	108	28	46	30	288	9302	2154	14254	16620
46	28	17532	29872	5846	51810	2346	938	286	22146	612	3346
3048	862	206	5668	3686	166	2092	214	4384	15762	14404	12586
390	64	3518	926	4362	318	176	1032	924	156	6956	2938
672	810	314	242	252	26	378	30	20	34	50	38

R5	No Drug										X4										No Drug										L83										AMD									
	34	14,140	46	21,656	30	3,988	30	49	42,206	45	44,589	32	26,471	55	6,849	38,144	7,336	31	1,217	1,128	70	1,159	97	29,118																										
	4,552	149	67,389	611	46	284	100	11,299	1,273	7,375	4,397	16,115	19,102	226	72	19,715	3,948	4,330	5,204	52	20	10,856	7,130																											
	67,828	14,982	3,788	50	658	231	53	177	24,739	612	2,997	416	68	12,020	4,384	37	242	35	974	32	33	45	9,695																											
	100	3,724	450	148	32	60	45	714	11,029	2,908	13,997	24,377	68	2,310	688	24	22	61	338	45	18	43	52																											
	35	83	18,099	39,257	5,413	97,733	9,254	5,249	470	19,175	513	3,264	107	4,209	8,451	38	40	29	23,846	20,393	39	126	31																											
	3,769	786	340	5,189	4,318	196	3,853	411	7,857	25,437	17,443	16,707	40	48	21	43	37	51	77	34	41	38																												
	391	98	4,449	4,357	6,090	1,886	263	3,890	2,089	475	8,475	4,107	26	38	28	9,102	105	3,635	38	9,577	1,698	43	38																											
	510	979	491	300	297	39	443	34	22	45	43	46	33	59	33	3,815	3,615	61,594	13,452	126,238	62	33	42	35																										
	L83																																																	
54		29	27	32	32	34	27	37	90	35	60	40	19,258	32	4,696	25,542	5,468	28	814	934	32	761	46	22,279																										
38		40	99	30	21	38	32	22	26	40	42	68	9,956	153	38	14,985	2,953	2,918	3,058	33	33	6,455	3,874																											
129		111	59	38	32	47	28	32	31	48	22	28	43	10,839	3,397	60	221	43	828	33	36	54	7,261																											
29		38	37	33	30	30	32	31	18	31	73	59	30	2,656	463	35	28	65	207	35	23	32	39																											
31		115	13,896	85	35	30	43	38	37	96	35	30	68	3,502	6,220	34	33	34	9,150	8,290	33	183	33																											
29		39	31	34	28	23	39	22	42	35	40	104	24	38	42	28	28	27	34	47	39	43																												
50		31	43	34	73	38	33	44	33	30	30	28	35	45	40	7,932	77	2,092	58	6,140	1,224	50	44																											
42		55	49	29	32	25	33	35	33	36	20	27	27	42	55	3,661	2,946	56,739	12,051	85,018	33	39	37	25																										
AMD																																																		
	39	12,186	32	7,226	36	3,664	34	40	33,192	44	30,510	35	32	28	40	32	39	142	42	37	40	33	37	31																										
	838	41	42,295	76	45	25	30	8,232	1,129	3,961	1,656	11,980	34	24	32	1,022	152	48	34	35	29	2,283	1,285																											
	51,146	8,580	2,049	53	493	219	26	361	21,767	479	925	597	39	3,384	538	49	38	34	95	28	34	40	35																											
	153	1,984	206	139	33	42	30	366	8,390	2,306	12,139	16,045	28	272	72	30	26	36	29	40	43	30	28																											
	38	3	16,980	27,832	5,043	48,452	2,219	784	310</																																									

Table 6

## *T20 Resistance Mutations*

<u>SDM Virus</u>	<u>T20 sens.<sup>a</sup></u>	<u>T20 FC<sup>b</sup></u>
HXB2 G I V		1.0
NL4-3 G I V	S	5.2
NL4-3 D I V	S	12.8
NL4-3 G I M	S	33.0
NL4-3 S I V	S	74.2
NL4-3 D I M	R	113.0
NL4-3 S I M	R	227.4
NL4-3 D T V	R	>281.8
JRCSF G I V		2.1
JRCSF D I V		104.0

<sup>a</sup> Rimsky et al., J. Virol. 72(2):986-993

<sup>b</sup> Fold change in IC50 (vs. HXB2) using PhenoSense HIV Entry Assay

### ***T20 Resistance Mutations***

<u>SDM Virus</u>	<u>T20 Sens.<sup>a</sup></u>	<u>T20 FC<sup>b</sup></u>
HXB2 G I V		1.0
NL4-3 G I V	S	5.2
NL4-3 D I V	S	12.8
NL4-3 G I M	S	33.0
NL4-3 S I V	S	74.2
NL4-3 D I M	R	113.0
NL4-3 S I M	R	227.4
NL4-3 D T V	R	>281.8
JRCSF G I V		2.1
JRCSF D I V		104.0

<sup>a</sup> Rimsky et al., J. Virol. 72(2):986-993

<sup>b</sup> Fold change in IC50 (vs. HXB2) using PhenoSense HIV Entry Assay

**Table 6**

#### **Entry Inhibitor Susceptibility: T-20 Fusion Inhibitor**

Drug resistance mutations were introduced into well-characterized X4 tropic (NL4-3) and R5 tropic (JRCSF) viruses. T20 susceptibility was measured using the virus entry assay. The fold change (FC) in T-20 susceptibility for each virus was determined by dividing the IC50 of the test virus by the IC50 of the HXB2 strain of HIV-1. T-20 sensitivity of similar mutant viruses has been reported in the scientific literature (Rimsky et al.).

# Identifying Entry Inhibitor Resistance Mutations

<u>Virus</u>	<u>AA Sequence<sup>a</sup></u>									
HXB2	Q	L	L	S	G	I	V	Q	Q	Q
1	Q	L	L	S	G	I	V	Q	Q	Q
2	Q	L	L	S		I		Q	Q	Q
3	Q	L	L	S	X	I	X	Q	Q	Q
4	Q	L	L	S	G	I	V	Q	Q	Q
5	Q	L	L	S	G	I	V	Q	Q	Q
6	Q	L	L	S	X	I	X	Q	Q	Q
7	Q	L	L	S	G	I	V	Q	Q	Q
8	Q	L	L	S	G	I	V	Q	Q	Q
9	Q	L	L	S	G	I	V	Q	Q	Q
10	Q	L	L	S		I	V	Q	Q	Q
12	Q	L	L	S	X	I	V	Q	Q	Q
13	Q	L	L	S	G	I	V	Q	Q	Q

X=G/S X=V/M

X=G/S X=V/M

X=G/D

<sup>a</sup> gp41 amino acid sequence positions 32 to 41

### Identifying Entry Inhibitor Resistance Mutations

<u>Virus</u>	<u>AA Sequence<sup>a</sup></u>		
HXB2	Q L L S G I V Q Q Q		
1	Q L L S G I V Q Q Q		
2	Q L L S S I M Q Q Q		
3	Q L L S X I X Q Q Q	X=G/S	X=V/M
4	Q L L S G I V Q Q Q		
5	Q L L S G I V Q Q Q		
6	Q L L S X I X Q Q Q	X=G/S	X=V/M
7	Q L L S G I V Q Q Q		
8	Q L L S G I V Q Q Q		
9	Q L L S G I V Q Q Q		
10	Q L L S D I V Q Q Q		
12	Q L L S X I V Q Q Q	X=G/D	
13	Q L L S G I V Q Q Q		

<sup>a</sup> gp41 amino acid sequence positions 32 to 41


**Table 7**

#### Identifying Envelope Protein Mutations.

Viruses with reduced (or increased) susceptibility to an entry inhibitor are identified using the virus entry assay. Mutations that may confer reduced (or increased) susceptibility to the entry inhibitor are identified by sequencing the envelope genes of the sensitive and resistant viruses. The deduced amino acid sequences of the sensitive and resistant viruses are compared to identify candidate drug resistance mutations. The ability of a specific mutation to confer altered drug susceptibility is confirmed or disproved by introducing the mutation into a drug sensitive virus and measuring the susceptibility of the mutant virus in the virus entry assay. In the example represented here, a short stretch of amino acid sequences within the first heptad repeat (HR-1) of the HIV-1 gp41 transmembrane envelope protein is aligned for viruses exhibiting different T-20 susceptibilities. Highlighted amino acids represent mutations known to confer reduced susceptibility to T-20.

Similar analyses can be used to identify envelope amino acid sequences that (a) alter/influence susceptibility to CCR5 or CXCR4 inhibitors, (b) specify X4, R5 and dual tropism, and © elicit neutralizing antibodies.

United States Patent & Trademark Office  
Office of Initial Patent Examination -- Scanning Division



Application deficiencies found during scanning:

☐ Page(s) \_\_\_\_\_ of \_\_\_\_\_ were not present  
for scanning. (Document title)

☐ Page(s) \_\_\_\_\_ of \_\_\_\_\_ were not present  
for scanning. (Document title)

They are 8 Drawings NOT 13 pages

☒ Scanned copy is best available. Drawings and Miscellaneous