

Table 1

# Cells

Cell	Receptor
5.25	CXCR4, CD4, CCR5 (not expressed well) BONZO
5.25.Luc4.M7	CD4, CCR5, BONZO
HOS.CD4.CCR5	CD4, CCR5
HOS.CD4.CXCR4	CD4, CXCR4
HOS.CD4	CD4, low level expression of CCR5 and CXCR4
HOS HT4 R5 GFP wt	CD4, CXCR4, CCR5
HOS.CD4.CCR5.GFP.M7#6*	CD4, CXCR4, CCR5
P4.CCR5	CD4, CXCR4, CCR5
U87.CD4	CD4
U87.CD4 R5	CD4, CCR5
U87.CD4 X4	CD4, CXCR4
MT2	CD4, CXCR4
MT4	CD4, CXCR4
PM1	CD4, CXCR4, CCR5
CEM NKr CCR5	CD4, CXCR4, CCR5

Table 2 Representative viruses and reagents

Viruses	Envelope <sup>a</sup>	Source
89.6, SF2	R5-X4/SI/B	ARRRP <sup>b</sup>
92BR014, 92US076	R5-X4/SI/B	ARRRP
JR-CSF, 91US005	R5/NSI/B	ARRRP
91US054	SI/B	ARRRP
NL43, MN, ELI	X4/B	ARRRP
92HT599	X4	ARRRP
92UG031	R5/NSI/A	ARRRP (IN-HOUSE)
92TH014, 92TH026	R5/NSI/B	ARRRP (IN-HOUSE)
92BR025, 93MW959	R5/SI/C	ARRRP (IN-HOUSE)
92UG035	R5/NSI/D	ARRRP (IN-HOUSE)
92TH022, 92TH023	R5/NSI/E	ARRRP (IN-HOUSE)
93BR020	R5-X4/SI/F	ARRRP (IN-HOUSE)
Antibodies	Epitope	SOURCE
Mabs 2F5, 1577	gp41 TM	ARRRP
Mabs IG1b12, 2G12, 17b, 48D	gp120 SU	ARRRP
Neutralization sera #2, HIV-IG	Polyclonal	ARRRP
Entry inhibitors	Target	Source
CD4-IG	gp120 SU	Genentech
CD4-IGG2	gp120 SU	Adarc
SCD4	Sigma	Progenics
T20 (DP178)	gp41 TM	Trimeris
Rantes, MIP1a/b	CCR5	SIGMA/ARRRP
SDF1a/b	CXCR4	SIGMA/ARRRP
AMD 3100	CXCR4	AnorMed
Dextran sulfate, Heparin	Non-specific	Sigma

<sup>a</sup>R5 (CCR5 co-receptor), X4 (CXCR4 co-receptor)

SI (syncytium inducing), NSI (non-syncytium inducing), A,B,C,D,E,F (envelope clade designation)

<sup>b</sup>AIDS Research and Reference Reagent Program



# Table 4 (Panel 1)

## Anti-HIV Drugs

Drug/Compound	Generic Name	Trademark	Manufacturer
<b>RT Inhibitors (NRTI, nucleotide analogs)</b>			
AZT, ZDV	Zidovudine	Retrovir	Glaxo/Wellcome
3TC	Lamivudine	Epivir	Glaxo/Wellcome
AZT + 3TC		Combivir	Glaxo/Wellcome
d4T	Stavudine	Zerit	Bristol-Myers/Squibb
ddI	Didanosine	Videx	Bristol-Myers/Squibb
ddC	Zalcitabine	Hivid	Hoffman La Roche
1592U89	Abacavir	Ziagen	Glaxo/Wellcome
AZT + 3TC + 1592U89		Trizivir	Glaxo/Wellcome
(-)FTC (5-fluoro-3TC; Corviracil)	Emtricitabine		Triangle Pharmaceuticals
(-)FTC + (+)FTC (50:50)	Racimir		QuadPharma
DAPD (DXG active)	Amdoxovir		Triangle Pharmaceuticals
F-ddA (2-fluoro-ddA)	Lodenosine		MedImmune Oncology (US Bioscience)
BCH-10652, dOTC (2-deoxy-3-oxa-4-thiocytidine)			BioChem Pharma, Inc.
D-d4FC			Triangle Pharmaceuticals (Schinazi)
<b>RT Inhibitors (NRTI, nucleotide analogs)</b>			
bis-POC PMPA (GS-4331)	Tenofovir		Gilead Sciences
bis-POM PMEa (GS-840)	Adefovir dipivoxil		Gilead Sciences
<b>RT Inhibitors (NNRTI, non-nucleosides)</b>			
BI-RG-587	Nevirapine	Viramune	Boehringer-Ingelheim (Roxanne)
BHAP PNU-90152T	Delavirdine	Rescriptor	Pharmacia & Upjohn
DMP 266 (L-743,726)	Efavirenz	Sustiva	Dupont Pharmaceuticals (Avid)
MKC442 (Coactinon)	Emivirine		Triangle/Mitsubishi Kasei
AG-1549 (S1153) (on hold)	Capravirine		Agouron Pharmaceuticals
PNU-142721			Pharmacia & Upjohn
DPC-961, -963, -083, -08?			DuPont Pharmaceuticals
SJ-3366	Also entry inhibitor?		Samjin Pharmaceuticals
BHAP PNU-87201	Atevirdine		Upjohn
GW420867X (quinoxaline)	(2 <sup>nd</sup> gen. HBY 097)		Glaxo/Wellcome (Hoechst Bayer)
TMC 120 (R147681)			Tibotec
TMC 125 (R165335)			Tibotec
R86183	tivirapine		Janssen Pharmaceuticals
Calanolide A			Sarawak Medichem Pharmaceuticals
<b>Protease Inhibitors (PRI)</b>			
Ro 31-8959	Saquinavir-(hgc) Saquinavir-(sgc)	Invirase Fortivase	Hoffman-La Roche
MK-639 (L-735,524)	Indinavir	Crixivar	Merck Research Laboratories
ABT-538 (A-84538)	Ritonavir	Norvir	Abbott Laboratories
AG1343	Nelfinavir	Viracept	Agouron Pharmaceuticals
141W94 (VX-478)	Amprenavir	Agenerase	Glaxo-Wellcome/Vertex
ABT-378/r	Lopinavir/ritonavir	Kaletra	Abbott Laboratories
BMS 232,632 (aza-peptide)			Bristol-Myers-Squibb
PNU-140690	Tipranavir		Pharmacia & Upjohn
DMP 450 (cyclic urea)	Mozenavir		Triangle/Avid (ph I/II)
TMC 126 (Erickson's compound)			Tibotec
GW433908 (VX-175)	amprenavir pro-drug		Glaxo/Wellcome/Vertex
L756,423 (on hold)			Merck
PD-178390 (dihydropyrone)			Parke Davis (Boehringer-Ingelheim)
? new candidate			Roche
DPC 681 and 684			DuPont Pharmaceuticals
AG-1776 (JE-2147 = KNI-764)			Agouron Pharmaceuticals
<b>Envelope/Receptor Inhibitors</b>			
T-20 (gp41)	Pentafuside		Trimeris Pharmaceuticals
T-1249 (gp41)			Trimeris Pharmaceuticals
D-peptide inhibitor (gp41) small mol.	SCH-C		Schering-Plough
AMD-3100 (CXCR4)	(bicyclam)		AnorMED
AMD-8664 (CXCR4)	(macrocyclam)		AnorMED
ALX40-4C (CXCR4)			U. PA
FP21399			Fuji Pharmaceuticals
PRO 542 (gp120)	CD4IgG2		Progenics Pharmaceuticals
PRO-140 (CCR5)	MAb CCR5		Progenics Pharmaceuticals
T-22 (CXCR4)	(peptide, 18-mer)		
Met-SDF-1 (CXCR4)			
TAK 779 (CCR5 antagonist)			Takeda
AOP-Rantes (CCR5)			Gryphon Sciences

# Table 4 (Panel 2)

Rantes 9-68 (CCR5)			
CCR5 antagonists	4-(piperidin-1-yl) butane class		Merck
$\alpha$ -Immunokine-NNS03 (CCR5, CXCR4)	$\alpha$ -cobratoxin		PhyloMed Corp.
<b>Integrase Inhibitors</b>			
AR-177	Zintevir		Aronex Pharmaceuticals
Diketo acids			Merck Research Laboratories
<b>Nucleocapsid Inhibitors</b>			
RB 2121	cyclic peptide p7 mimic		(see PNAS 96:4886-4891 (1999))
CI-1012			Achelion Pharmaceuticals
<b>RNAse H Inhibitor</b>			
SP1093V (BBNH Fe+3 derivative)			(Parniak)

FDA approved drugs are shown in boldface, red = discontinued development, blue = not sure about development status

**Table 4 (Panel 1)**

<b>Generic Name (abbreviation)</b>	<b>Brand Name</b>	<b>Firm</b>	<b>FDA Approval Date</b>
zidovudine, AZT	Retrovir	Glaxo Wellcome	March 87
didanosine, ddI	Videx	Bristol Myers-Squibb	October 91
zalcitabine, ddC	Hivid	Hoffman-La Roche	June 92
stavudine, d4T	Zerit	Bristol Myers-Squibb	June 94
lamivudine, 3TC	Epivir	Glaxo Wellcome	November 95
saquinavir, SQV, hgc	Invirase	Hoffman-La Roche	December 95
saquinavir, SQV, sgc	Fortovase	Hoffman-La Roche	November 97
ritonavir, RTV	Norvir	Abbott Laboratories	March 96
indinavir, IDV	Crixivan	Merck & Co., Inc.	March 96
nevirapine, NVP	Viramune	Boehringer Ingelheim	June 96
nelfinavir, NFV	Viracept	Agouron Pharmaceuticals	March 97
delavirdine, DLV	Rescriptor	Pharmacia & Upjohn	April 97
ZDV+3TC	Combivir	Glaxo Wellcome	September 97
efavirenz, EFV	Sustiva	DuPont Pharmaceuticals	September 98
abacavir, ABC	Ziagen	Glaxo Wellcome	February 99
amprenavir	Agenerase	Glaxo Wellcome	April 99
lopinavir/ritonavir	Kaletra	Abbott	September 2000
ZDV+3TC+ABC	Trizivir	GlaxoSmithKline	November 2000

5.

drug 1: L83  
drug 2: AMD

	1			2			3			4			5			6		
No Drug RLU L83 RLU AMD RLU	R5 34	X4 26,471 19,258 32	R5:X4 0.00	R5 14,140 54 12,186	X4 55	R5:X4 257.09	R5 46	X4 6,849 4,696 40	R5:X4 0.01	R5 21,656 38,144 27 7,226	X4 14 25,542 32	R5:X4 0.57	R5 30	X4 7,336 5,468 39	R5:X4 0.00	R5 3,988 34 3,664	X4 31	R5:X4 128.65
No Drug RLU L83 RLU AMD RLU	R5 4,552 28 838	X4 19,102 9,956 34	R5:X4 0.24	R5 149 40 41	X4 226 153 24	R5:X4 0.66	R5 67,389 99 42,295	X4 72	R5:X4 935.96	R5 611 30 76	X4 19,715 14,985 1,022	R5:X4 0.03	R5 46	X4 3,948 2,953 152	R5:X4 0.01	R5 284 38 25	X4 4,330 2,918 48	R5:X4 0.07
No Drug RLU L83 RLU AMD RLU	R5 67,828 129 51,146	X4 68 997.47	R5:X4	R5 14,982 111 8,580	X4 12,020 10,839 3,384	R5:X4 1.25	R5 3,768 59 2,049	X4 4,384 3,397 538	R5:X4 0.86	R5 50 37	X4 37	R5:X4	R5 658 32 493	X4 242 221 38	R5:X4 2.72	R5 231 47 219	X4 35	R5:X4 6.60
No Drug RLU L83 RLU AMD RLU	R5 100	X4 39	R5:X4	R5 3,724 38 1,984	X4 2,310 2,656 272	R5:X4 1.61	R5 450 37 206	X4 668 483 72	R5:X4 0.67	R5 148 33 139	X4 24	R5:X4 6.17	R5 32	X4 22	R5:X4	R5 60	X4 61	R5:X4
No Drug RLU L83 RLU AMD RLU	R5 35	X4 107 68 34	R5:X4 0.33	R5 83	X4 4,209 3,502 88	R5:X4 0.02	R5 18,099 13,896 16,980	X4 8,451 6,220 8,711	R5:X4 2.14	R5 39,257 85 27,832	X4 38 1033.08	R5:X4	R5 5,413 35 5,043	X4 40 135.33	R5:X4	R5 97,733 30 48,452	X4 29 3370.10	R5:X4
No Drug RLU L83 RLU AMD RLU	R5 3,769 29 2,773	X4 40 94.23	R5:X4	R5 786 39 797	X4 48 16.38	R5:X4	R5 340 31 282	X4 21 16.19	R5:X4	R5 5,189 34 5,356	X4 43 120.67	R5:X4	R5 4,318 28 3,461	X4 37 118.70	R5:X4	R5 196 23 184	X4 51	R5:X4 3.84
No Drug RLU L83 RLU AMD RLU	R5 391 50 340	X4 26 15.04	R5:X4	R5 98 38	X4 38	R5:X4	R5 4,449 43 3,752	X4 28 7,932 906	R5:X4 158.89	R5 4,357 34 906	X4 9,102 7,932 36	R5:X4 0.48	R5 6,090 73 4,473	X4 105 77 51	R5:X4 58.00	R5 1,886 38 279	X4 3,635 2,092 27	R5:X4 0.52
No Drug RLU L83 RLU AMD RLU	R5 510 42 564	X4 33 15.45	R5:X4	R5 979 55 813	X4 59 16.59	R5:X4	R5 491 49 376	X4 33 14.88	R5:X4	R5 300 29 268	X4 3,815 3,661 36	R5:X4 0.08	R5 297 32 256	X4 3,615 2,946 30	R5:X4 0.08	R5 39 56,739 40	X4 61,594 56,739 40	R5:X4 0.00

[illegible]

1

	3		
	R5 cells	X4 cells	R5:X4
no drug			
R5 inhibitor			
X4 inhibitor			
%inhib by R5 inhibitor			
%inhib by X4 inhibitor			

	21		
	R5	X4	R5:X4
no drug			
R5 inhibitor			
X4 inhibitor			
%inhib by R5 inhibitor			
%inhib by X4 inhibitor			

[illegible]

## DUAL or MIXED

DEAD

	76		
	R5	X4	R5:X4
no drug	14,982	12,020	1
R5 inhibitor	111	10,839	
X4 inhibitor	8,580	3,384	
%inhib by R5 inhibitor	99	10	
%inhib by X4 inhibitor	43	72	

	R5	X4	R5:X4
no drug	43	42	
R5 inhibitor			
X4 inhibitor			
%inhib by R5 inhibitor			
%inhib by X4 inhibitor			

[illegible]



Table

Plate Repeat End time Start temp. End temp. BarCode  
1 1 6:39:36 PM 21.6 21.7 N/A

0.5 CPS (CPS)  
18

26648	54	6970	37406	7158	38	1386	930	100	1184	112	30488
19248	262	72	18972	3478	3946	5010	46	14	11004	7474	60
106	9038	4002	32	238	26	976	34	38	48	9736	32
38	2000	528	16	16	36	298	52	14	54	52	38
122	2984	7264	40	38	20	23344	15340	44	156	40	32
36	52	28	40	28	40	62	32	38	22	44	28
22	32	28	9858	68	3802	46	11470	1958	42	42	24
30	64	36	3846	3390	57858	12620	126186	68	34	38	36

Plate Repeat End time Start temp. End temp. BarCode  
2 1 6:41:51 PM 21.5 21.6 N/A

0.5 CPS (CPS)  
34

26294	56	6728	38882	7514	24	1048	1326	40	1134	82	27748
18956	190	72	20458	4418	4714	5398	58	26	10708	6786	40
30	15002	4766	42	246	44	972	30	28	42	9654	30
40	2620	808	32	28	86	378	38	22	32	52	30
92	5434	9638	36	42	38	24348	25446	34	96	22	24
44	44	14	46	46	62	92	36	44	54	46	48
30	44	28	8346	142	3468	30	7684	1438	44	34	42
36	54	30	3784	3840	65330	14284	130290	56	32	46	34

Plate Repeat End time Start temp. End temp. BarCode  
3 1 6:44:06 PM 21.6 21.6 N/A

0.5 CPS (CPS)  
16

18590	18	4306	23902	5386	30	924	894	32	660	48	23382
8698	148	34	14088	2880	3142	3160	40	46	7616	2842	32
30	10252	3542	38	172	34	842	30	32	44	7616	34
28	2396	370	28	26	66	172	32	22	28	36	26
64	3784	4822	32	28	24	6188	5702	34	120	34	24
16	28	62	26	28	30	38	38	38	44	62	34
38	52	40	7020	68	1798	62	7324	1076	50	42	32
32	50	70	3824	3138	53670	13088	104608	32	38	50	32

Plate Repeat End time Start temp. End temp. BarCode  
4 1 6:46:20 PM 21.7 21.8 N/A

0.5 CPS (CPS)  
24

19926	46	5086	27182	5550	26	704	974	32	862	44	21176
11214	158	42	15882	3026	2694	2956	26	20	5294	4906	44
56	11426	3252	82	270	52	814	36	40	64	6906	32
32	2916	556	42	30	64	242	38	24	36	42	40
72	3220	7618	36	38	44	12112	10878	32	246	32	42
32	48	22	30	28	24	136	30	56	34	24	24
32	38	40	8844	86	2386	54	4956	1372	50	46	28
22	34	40	3498	2754	59808	11014	65428	34	40	24	18

Plate Repeat End time Start temp. End temp. BarCode  
5 1 6:48:35 PM 21.9 21.9 N/A

0.5 CPS (CPS)  
38

38	28	40	30	38	234	46	30	46	30	34	22
26	20	32	976	190	58	46	24	28	2298	1290	26
42	3260	420	50	36	34	102	26	34	24	40	24
28	378	88	32	34	32	40	42	54	26	30	48
36	62	11690	38	44	32	42	42	26	36	46	42
46	22	44	38	58	32	38	30	32	36	38	12
22	38	52	40	42	26	34	50	200	36	34	36
36	80	52	40	36	40	52	54	22	40	42	36

Plate Repeat End time Start temp. End temp. BarCode  
6 1 6:50:50 PM 21.8 21.9 N/A

0.5 CPS (CPS)  
30

26	28	40	34	40	50	38	44	34	36	40	40
42	28	32	1068	114	38	22	46	30	2268	1280	40
36	3508	656	48	40	34	88	30	34	56	30	36
28	166	56	28	18	40	18	38	32	34	26	32
32	114	5732	32	42	34	38	34	50	34	38	52
36	22	36	40	30	20	30	26	42	40	26	36
46	24	26	32	60	28	32	40	120	36	20	36
32	58	58	32	24	40	50	36	12	38	40	18

Talla

Plate Repeat End time Start temp. End temp. BarCode  
1 1 ##### 21.7 21.8 N/A

0.5 CPS (CPS)  
32

24	14218	56	21586	28	4034	36	52	42088	40	43332	38
4780	174	71276	730	48	326	28	12022	1288	7198	4314	17856
69452	15306	4008	44	792	180	50	182	25292	718	3572	424
72	3922	546	128	44	66	30	900	11984	3194	16106	23794
30	88	18820	40438	4882	105946	8466	4934	470	20386	548	3376
3420	934	344	6268	5012	170	3546	504	8164	22214	15146	19592
344	82	5344	5330	6710	1880	338	4286	2112	466	7752	4078
580	1018	516	318	302	40	414	32	20	52	48	42

Plate Repeat End time Start temp. End temp. BarCode  
2 1 ##### 21.9 21.9 N/A

0.5 CPS (CPS)  
22

44	14062	36	21726	32	3942	24	46	42324	50	45846	26
4324	124	63502	492	44	242	172	10576	1258	7552	4480	14374
66204	14658	3568	56	524	282	56	172	24186	506	2422	408
128	3526	354	168	20	54	60	528	10074	2622	11888	24960
40	78	17378	38076	5944	89520	10042	5564	470	17964	478	3152
4118	638	336	4110	3624	222	4160	318	7550	28660	19740	13822
438	114	3554	3384	5470	1892	188	3494	2066	484	9198	4136
440	940	466	282	292	38	472	36	24	38	38	50

Plate Repeat End time Start temp. End temp. BarCode  
3 1 ##### 21.9 22 N/A

0.5 CPS (CPS)  
36

36	38	28	28	38	28	36	42	114	40	60	46
32	34	80	26	16	36	30	18	26	30	48	80
126	106	34	42	30	46	22	26	32	58	8	28
26	46	42	38	34	20	18	34	24	36	100	60
28	146	15756	76	32	26	42	52	36	106	34	24
34	50	20	24	32	18	38	18	36	38	50	58
38	38	28	36	80	50	38	44	34	38	30	28
40	30	38	30	24	22	38	36	34	36	22	26

Plate Repeat End time Start temp. End temp. BarCode  
4 1 ##### 21.9 22 N/A

0.5 CPS (CPS)  
16

24	70	30	26	26	40	18	32	66	30	60	34
24	46	118	34	26	40	34	26	26	50	36	56
132	116	84	34	34	48	34	38	30	38	36	28
32	30	32	28	26	40	46	28	12	26	46	58
34	84	12036	94	38	34	44	24	38	86	36	36
24	28	42	44	24	28	40	26	48	32	30	150
62	24	58	32	66	26	28	44	32	22	30	28
44	80	60	28	40	28	28	34	32	36	18	28

Plate Repeat End time Start temp. End temp. BarCode  
5 1 ##### 22 22 N/A

0.5 CPS (CPS)  
34

38	11132	26	6696	34	3960	44	40	29548	40	30628	34
762	48	39888	88	40	30	36	8416	1262	4096	1598	12372
52058	8050	2524	60	450	168	32	340	22410	534	1046	512
106	2080	172	170	38	38	30	444	7478	2458	10024	15470
30	46	16428	25792	4240	45094	2092	630	334	17130	388	2730
2498	732	358	5044	3236	202	2124	292	4806	13736	13012	12386
290	30	3986	886	4584	240	110	940	1286	162	6768	2204
456	816	438	294	260	44	524	40	34	22	38	20

Plate Repeat End time Start temp. End temp. BarCode  
6 1 ##### 22 22.1 N/A

0.5 CPS (CPS)  
32

40	13240	38	7756	38	3368	24	40	36836	48	30392	36
914	34	44702	64	50	20	24	8048	996	3826	1714	11588
50234	9110	1574	46	536	270	20	382	21124	424	804	682
200	1888	240	108	28	46	30	288	9302	2154	14254	16620
46	28	17532	29872	5846	51810	2346	938	286	22146	612	3346
3048	862	206	5668	3686	166	2092	214	4384	15762	14404	12586
390	64	3518	926	4362	318	176	1032	924	156	6956	2938
672	810	314	242	252	26	378	30	20	34	50	38

[illegible]

Table 6

# T20 Resistance Mutations

<u>SDM virus</u>	<u>T20 Sens.<sup>a</sup></u>	<u>T20 FC<sup>b</sup></u>
HXB2 G I V		1.0
NL4-3 G I V	S	5.2
NL4-3 D I V	S	12.8
NL4-3 G I M	S	33.0
NL4-3 S I V	S	74.2
NL4-3 D I M	R	113.0
NL4-3 S I M	R	227.4
NL4-3 D T V	R	>281.8
JRCSF G I V		2.1
JRCSF D I V		104.0

<sup>a</sup> Rimsy et al., J. Virol. 72(2):986-993

<sup>b</sup> Fold change in IC50 (vs. HXB2) using Phenosense HIV Entry Assay

### ***T20 Resistance Mutations***

<u>SDM Virus</u>	<u>T20 Sens.<sup>a</sup></u>	<u>T20 FC<sup>b</sup></u>
HXB2 G I V		1.0
NL4-3 G I V	S	5.2
NL4-3 D I V	S	12.8
NL4-3 G I M	S	33.0
NL4-3 S I V	S	74.2
NL4-3 D I M	R	113.0
NL4-3 S I M	R	227.4
NL4-3 D T V	R	>281.8
JRCSF G I V		2.1
JRCSF D I V		104.0

<sup>a</sup> Rimsky et al., J. Virol. 72(2):986-993

<sup>b</sup> Fold change in IC50 (vs. HXB2) using PhenoSense HIV Entry Assay

**Table 6**

#### **Entry Inhibitor Susceptibility: T-20 Fusion Inhibitor**

Drug resistance mutations were introduced into well-characterized X4 tropic (NL4-3) and R5 tropic (JRCSF) viruses. T20 susceptibility was measured using the virus entry assay. The fold change (FC) in T-20 susceptibility for each virus was determined by dividing the IC50 of the test virus by the IC50 of the HXB2 strain of HIV-1. T-20 sensitivity of similar mutant viruses has been reported in the scientific literature (Rimsky et al.,).

# Identifying Entry Inhibitor Resistance Mutations

<u>Virus</u>	<u>AA Sequence<sup>a</sup></u>									
HXB2	Q	L	L	S	G	I	V	Q	Q	Q
1	Q	L	L	S	G	I	V	Q	Q	Q
2	Q	L	L	S	I			Q	Q	Q
3	Q	L	L	S	X	I	X	Q	Q	Q
4	Q	L	L	S	G	I	V	Q	Q	Q
5	Q	L	L	S	G	I	V	Q	Q	Q
6	Q	L	L	S	X	I	X	Q	Q	Q
7	Q	L	L	S	G	I	V	Q	Q	Q
8	Q	L	L	S	G	I	V	Q	Q	Q
9	Q	L	L	S	G	I	V	Q	Q	Q
10	Q	L	L	S	I	V		Q	Q	Q
12	Q	L	L	S	X	I	V	Q	Q	Q
13	Q	L	L	S	G	I	V	Q	Q	Q

X=G/S X=V/M

X=G/S X=V/M

X=G/D

<sup>a</sup> gp41 amino acid sequence positions 32 to 41

### Identifying Entry Inhibitor Resistance Mutations

<u>Virus</u>	<u>AA Sequence<sup>a</sup></u>	
HXB2	Q L L S G I V Q Q Q	
1	Q L L S G I V Q Q Q	
2	Q L L S S I M Q Q Q	
3	Q L L S X I X Q Q Q	X=G/S X=V/M
4	Q L L S G I V Q Q Q	
5	Q L L S G I V Q Q Q	
6	Q L L S X I X Q Q Q	X=G/S X=V/M
7	Q L L S G I V Q Q Q	
8	Q L L S G I V Q Q Q	
9	Q L L S G I V Q Q Q	
10	Q L L S D I V Q Q Q	
12	Q L L S X I V Q Q Q	X=G/D
13	Q L L S G I V Q Q Q	

<sup>a</sup> gp41 amino acid sequence positions 32 to 41

**Table 7**

### Identifying Envelope Protein Mutations.

Viruses with reduced (or increased) susceptibility to an entry inhibitor are identified using the virus entry assay. Mutations that may confer reduced (or increased) susceptibility to the entry inhibitor are identified by sequencing the envelope genes of the sensitive and resistant viruses. The deduced amino acid sequences of the sensitive and resistant viruses are compared to identify candidate drug resistance mutations. The ability of a specific mutation to confer altered drug susceptibility is confirmed or disproved by introducing the mutation into a drug sensitive virus and measuring the susceptibility of the mutant virus in the virus entry assay. In the example represented here, a short stretch of amino acid sequences within the first heptad repeat (HR-1) of the HIV-1 gp41 transmembrane envelope protein is aligned for viruses exhibiting different T-20 susceptibilities. Highlighted amino acids represent mutations known to confer reduced susceptibility to T-20.

Similar analyses can be used to identify envelope amino acid sequences that (a) alter/influence susceptibility to CCR5 or CXCR4 inhibitors, (b) specify X4, R5 and dual tropism, and © elicit neutralizing antibodies.